

Ocular Sagittal Height as a Better Predictor of Scleral Lens Sagittal Depth than Keratometry in a Case of Post-LASIK Ectasia

Claire Chu OD, Mari Fujimoto OD FAAO Dipl AAO, Patrick Caroline FAAO, Matthew Lampa OD FAAO, Mary Avisado OD FAAO, Randy Kojima FAAO
Pacific University College of Optometry, Forest Grove, Oregon

Introduction

Post-LASIK ectasia can cause irregular astigmatism that decreases spectacle best-corrected visual acuity (BCVA). Rigid contact lenses, including gas-permeable and scleral lenses (SL), often provide better vision for patients.¹

In a diagnostic fit, the initial SL is often selected based on keratometry (K) values, with steeper Ks indicating a deeper sagittal depth. However, K values alone may not allow for accurate diagnostic SL selection in an ectatic cornea.

Case Report

A 49-year-old systemically healthy White male was referred for contact lens fit for post-LASIK ectasia OD>OS. He had LASIK OU about 25 years ago and started habitually wearing spectacle correction one month ago. He has no prior history of contact lens wear.

Incoming spectacle VA was 20/400 OD PH NI, 20/25 OS PH 20/20. Slit lamp exam revealed clear central corneas with LASIK flap scars OU. Slit lamp exam, pachymetry, and AS-OCT demonstrated central corneal steepening and thinning OD.




Figure 1. Topographies OU
OD Ks 66.34D@173/72.60D@083
 Corneal astigmatism 6.26D
 Spectacle Rx +0.75-4.25x073
OS Ks 42.55D@175/45.54D@085
 Corneal astigmatism 2.99D
 Spectacle Rx -0.50-2.50x097

Post-LASIK Ectasia

- Loss of biomechanical strength leading to progressing myopia or astigmatism, corneal steepening, and corneal thinning¹
- Risk factors include abnormal topography, multiple procedures, high myopia, thin corneas, and thin residual stromal bed²⁻⁴
- Surgical management includes intracorneal ring segments, corneal cross-linking, and PKP or DALK^{1,3}
- Tomography may provide additional data for ectasia risk assessment, especially in conjunction with AI.^{3,4}

Contact Lens Management

AS-OCT demonstrated ocular sag at a 15.0mm chord was 3900µm OD. Despite presenting with steep corneal curvatures, a 16.5mm diagnostic scleral lens with a 4200 sag provided appropriate apical clearance of 300µm. Over-refraction of pl-2.50x150 improved vision to 20/20.

Patient elected to continue spectacle correction OS.




Figure 4. Ocular sag at a 15.0mm chord was 3900µm OD, 3260µm OS




Figure 2. Pachymetry demonstrates corneal thinning OD. Central corneal thickness 313µm OD, 473µm OS. Minimum thickness 250µm OD, 444µm OS.




Figure 3. AS-OCT shows thinning and hyperreflectivity of the central stroma OD.

Final Scleral Lens

OD: SL/8.04/-3.50-2.25x001/16.5/4200 sag/100µm toric haptic/Dk 100
VA: 20/20




Figure 5. AS-OCT demonstrating scleral lens clearance OD post-settling.




Figure 6. Final scleral lens OD.

Conclusion

In this case of post-LASIK ectasia, the steep K of 72.60D OD suggests that the patient may need a SL of exceptional depth to vault the cornea. However, ocular sag measurement revealed that the patient had sagittal height only slightly greater than the normal eye average of about 3700µm⁵, so a SL with average sag was appropriate for the patient. In SL fits for post-surgical ectasia, ocular sag measurements may more accurately predict the appropriate SL sag than Ks.

References

- Wolle MA, Randleman JB, Woodward MA. Complications of Refractive Surgery: Ectasia After Refractive Surgery. *Int Ophthalmol Clin.* 2016;56(2):127-39.
- Bohac M, Koncarevic M, Pasalic A, Biscevic A, Merlak M, Gabric N, Patel S. Incidence and Clinical Characteristics of Post LASIK Ectasia: A Review of over 30,000 LASIK Cases. *Semin Ophthalmol.* 2018;33(7-8):869-877.
- Zhao L, Yin Y, Hu T, Du K, Lu Y, Fu Q, Zhang Y, Wu X, Li Y, Wen D. Comprehensive management of post-LASIK ectasia: From prevention to treatment. *Acta Ophthalmol.* 2023;101(5):485-503.
- El-Naggar MT, Elkitkat RS, Ziada HE, Esporcatte LPG, Ambrósio R Jr. Assessment of Preoperative Risk Factors for Post-LASIK Ectasia Development. *Clin Ophthalmol.* 2023 4;17:3705-3715.
- Rojas Viñuela J, Piñero DP, Martínez MB. Comparing sagittal heights calculated using corneal parameters and those measured with profilometry. *Cont Lens Anterior Eye.* 2023;46(2):101747.