



Ampleye Standard 16.5 mm Diagnostic Set - Item# R1602

9 Ampleye Dx Lenses with 16.5mm diameter and toric haptic. *INCLUDES:* Sagittal depths of 3800, 4000, 4200, 4400, 4600, 4800, 5000, 5200, & 5600.

Ampleye 15.0/15.5 Diagnostic Set - Item# R1603

4 Ampleye Dx Lenses with 15.0mm diameter and 5 Ampleye Dx Lenses with 15.5mm diameter and toric haptic. *INCLUDES:* Sagittal depths of 3400, 3800, 4200, 4600 in both diameters and 4800.

Ampleye Standard 15.5 mm Diagnostic Set - Item# R1603B

9 Ampleye Dx Lenses with 15.5mm diameter and toric haptic. *INCLUDES:* Sagittal depths of 3200, 3400, 3600, 3800, 4000, 4200, 4400, 4600, & 4800.

Ampleye 15.5/16.5 Combo Diagnostic Set - Item# R1605

18 Ampleye Dx Lenses with 15.5/16.5 mm diameter and toric haptic. *INCLUDES:* 1 lens of 3200, 3400, 3600, 3800, 4000, 4200, 4400, 4600, 4800 sags in 15.5 diameter, as well as one each of sag depths 3800, 4000, 4200, 4400, 4600, 4800, 5000, 5200, & 5600 in 16.5 mm.

The above listed sets ship in storage case with fitting instructions and accessories - large DMV, small DMV, hand-held cobalt blue light, wratten filter, & starter supply of non-preserved saline.

Ampleye Diagnostic Set Extensions - Item# R1602EXT or R1603EXT

Sold as add-ons to Ampleye R1602 or R1603 9-lens sets. **R1602EXT** includes duplicate 16.5mm lenses with 3800, 4000, 4200, 4400, 4800 as well as a 5400sag, all with toric haptic, and spherical lenses of 4000, 4200 & 4400sag. **R1603EXT** includes 9 Ampleye Dx Lenses with 15.5 mm diameter and toric haptic in additional sag depths of 3200, 3400, 3600, 3800, 4000, 4200, 4400, 4600. Lenses will ship with an 18-lens foam insert that fits the original storage case.

Ampleye Individual Trial Lens - Item# R1602IRL

Available to order by diameter & sag for existing set-holders.



**Ampleye Patient Starter Kit
Item# R1130**

INCLUDES: deep well lens case, large DMV for lens insertion, traditional DMV for lens removal, pocket mirror, starter supply of non-preserved saline, & instructions in clear carrying case.

Progent by Menicon

For biweekly cleaning of scleral and corneal GP lenses. Cleans, loosens and removes surface deposits after a short soak without manual rubbing.

INCLUDES: Progent & lens case.



Item# R1730 - Single Treatment

Item#R1732 - 7 Treatment Pack



**Progent by Menicon Scleral Lens Case
Item# R1746**

For Progent treatment of large diameter lenses (11-23mm).

**Menicon LacriPure Saline Solution
Item# R1127**

98 Count Box - 5ml vials. FDA approved non-preserved saline solution indicated for rinsing soft, hybrid, or GP lenses and for scleral lens insertion.



**DMV Traditional RGP Lens Removers
Item# R3100**

Small suction cup for removing gas permeable lenses from the eye. Instructions included.



DMV Scleral Cup - Item# R3130

Large diameter GP lens insertion tool.



**EZI Scleral Lens Applicator
Item# R3135**

Scleral lens insertion ring allows for one finger application, lens self-positioning, and less air entrapment.

